

# Presentation 6: PFGE Troubleshooting

Everything you wanted to know about  
PFGE, but were afraid to ask...

**OR**

Yes, in fact, we've done that  
experiment!

# Trouble shooting philosophy

- ❑ Can or do we **test** every issue / variable?
  - Most
- ❑ Can or do we **solve** every issue / variable?
  - No, only consistent issues, but there are exceptions...
- ❑ Our role at PFGE Reference lab at the CDC
  - Test the big issues
  - Make recommendations

## Your expectations of us:

- ✓ get trained
- ✓ receive prompt answers to questions
- ✓ get feedback on gels
- ✓ receive occasional trouble shooting assistance
- ✓ read about reviews and reports on reagent options

## Our expectations of labs:

- ✓ **follow protocols closely**
- ✓ perform trouble shooting
- ✓ provide feedback
- ✓ **follow protocols closely**

# Trouble shooting approach

- ❑ **What has changed since the last “good” gel?**
  - Equipment, reagent, person, isolate / strain, etc...
  - Change in lab procedures / new SOPs, new vendor, contamination / decontamination event, etc...
- ❑ **Examine the gel closely**
  - Ask if the problem is:
    - in all lanes, only the isolates, only the standards, only one enzyme?
    - apparent with multiple lab personnel?
    - apparent with multiple enzymes or multiple organisms?
    - associated with a specific mapper?
- ❑ **Maintain notebook of troubleshooting findings, gels**
  - ➔ **What is in common and what has changed?**

# Trouble shooting approach, continued

## ❑ Reagents: garbage in → garbage out

- Always track vendors or individuals (if in-house), dates, lots, etc...
- Make and store properly – room temperature, 4°C / -20°C, protect from light, dilute in TE or water, consider making aliquots, etc...
- If it looks suspicious → do not use it!
  - cloudy, stringy, precipitates, etc...
  - do not be superstitious
- Adverse affects may be gradual or sudden, may be pronounced in some organisms and not others or H9812




## ❑ Difficulty repeating these effects and “bad gels”

- Complex experiment design, interaction of variables, protocol interpretation, unique lab environment, too many variables!!!
- All organisms, all strains, all enzymes not affected equally


## Example of bad cell treatment – plates grown for too long

18 hrs



yes

48 hrs




no

- ❑ Plates should be grown for **14 – 18 hours**
- ❑ **Fresh** cell growth
- ❑ Do not keep plates at 4°C, grow on bench, let sit out at room temperature, etc...
- ❑ Treat cells carefully
- ❑ Not all strains / serotypes similarly affected

# Example of bad cell treatment – plates stored improperly

*E. coli* O157 cultures stored at 4°C for 3 days before making cell suspensions




no


- ❑ Plates should be grown for **14 – 18 hours**
- ❑ **Fresh** cell growth
- ❑ Do not keep plates at 4°C, grow on bench, let sit out at room temperature, etc...
- ❑ Treat cells carefully
- ❑ Not all strains / serotypes similarly affected

# Thiourea-dependent or “Untypeable” Strains

No thiourea



50 $\mu$ M thiourea




- Typical smear – DNA degraded in bottom third of gel, no DNA in wells
- Add thiourea to the running buffer
  - 1g thiourea in 100ml dH<sub>2</sub>O (10mg/ml)
  - wrap bottle in foil and store in dark
  - add 837 $\mu$ L to 2.2L 0.5X TBE at start of run (50mM)
- Do not use thiourea routinely or add directly to plugs because it is toxic
- Thiourea does not improve poor quality plugs or poor patterns
- Include a positive control on gel
  - digested plug slice of “untypeable” strain that “typed” previously when thiourea was added to buffer

# Example of bad plug – thiourea will not help

## □ Are these thiourea-dependent isolates?


- Is there DNA in the plug slice?
- Is the smear in the whole lane or just towards the bottom?
- Are there any bands present?

### No thiourea in running buffer



thiourea-dependent


### 50 $\mu$ M thiourea in running buffer



bad plug



# Incorrect stopping point – plugs stored before lysis




**What to consider** → lysis buffer, cell suspension buffer, reagent contamination, reagent storage, cell “happiness,” alteration of standardized protocol

# Example of contaminated mapper


## ❑ Contamination in the mapper (fungal, bacterial)

- Smearing throughout the lane (usually)
- Smearing in isolates and controls (usually)

**Contaminated mapper**



**Same mapper after being decontaminated**




# Example of Contaminated Mapper

## ❑ Solution

- Drain 0.5X TBE / water from mappers / lines after **each use**
- Once every 3 months (or as needed)
  - circulate ~2 L 5% bleach 20 – 30 minutes, drain
  - circulate ~2 L dH<sub>2</sub>O 20 – 30 minutes, drain, repeat at least 3X
  - leave chiller OFF

Pseudomonas species isolated  
from CHEF Mapper unit at CDC




# Example of Contaminated Plugs

## ❑ Bad *Xba*I / H buffer?

- Unlikely → completely degraded DNA, no ghost bands


## ❑ Solution

- Keep plugs in sterile TE at 4°C
- Keep track of dates – identify and discard “bad” plugs
- Make a new batch of H9812 every 1 – 3 months, depending on testing volume, storage space
- Always use previously tested lot of H9812
- Make H9812 control plug
  - make new H9812 plug same day as unknowns
  - controls for plug making that day




# Troubleshooting missing bands

- **H9812 is the universal standard used with all organisms**
  - chosen due to optimal number and distribution of bands
- **Always test “new” lot of H9812 plugs with “old” lot to confirm:**
  - PFGE pattern is same, including band intensity
  - concentration is similar
  - the new standards produce good quality PFGE pattern
  - if pattern is different, streak out new plate from frozen stock
- **Pre-tested H9812 standard plugs can be a positive control and focus troubleshooting efforts**
  - problem observed in sample lanes but **not** in standards, likely due to sample associated affect (i.e. sample plug preparation)
  - problem observed in both sample **and** standard lanes, likely due to step involving sample and standard plugs (i.e. enzyme digests)



# Trouble shooting no bands – no digestion


- Some strains do not restrict with certain enzymes – no site or sites are “hidden”
- Appears as a large band at the top of the gel with no DNA smearing in the lane and DNA still left in the well
- It's not you – it's the strain (and enzyme)!
  - *E. coli*, *Salmonella*, *Listeria*, *Vibrio spp*
  - some serotypes more often than others
  - cannot predict



# Trouble shooting no bands – no digestion

- Other slices digested with master mix are a positive control for enzyme working
- Valid pattern for that strain and can be uploaded
  - not the same as untypeable!
- **Solution → digest with other enzymes (secondary +/- tertiary) for that organism**

*Listeria, Apal*





# Troubleshooting wavy bands

- **Possible causes**



- rusty, nicked, damaged razor or scalpel
- damaged plug slices, debris in gel and/or plug

- **Possible solutions**

- use new razor or scalpel
- make new plugs
- completely melt plug agarose and mix completely with cell suspension
- completely melt gel agarose, check for cell debris
- **do not** use alcohol wipes between plug slices

Spatula wiped with alcohol between handling digested slices


No alcohol





# “Ghost” or “Shadow” Bands

- **Due to incomplete digestion or star activity**
- **May be the result of:**
  - Poor plug quality
    - proteinase K not washed out of plug
    - enzyme inhibitors not washed out of plug
    - cell concentration too high (DNA and debris)
  - Poor enzyme and/or buffer quality
    - bad lot, change in manufacturing process
    - expired or vial opened frequently
  - Enzyme digestion not optimal
    - old/bad BSA or BSA not included in master mix
    - not enough units of enzyme
    - too many units of enzyme (star activity)
    - incubation time too short




- incubation time too long (star activity)
- incorrect incubation temperature
- incorrect buffer

# Troubleshooting “ghost” bands

- **Possible solutions**

- decrease cell suspension concentration
- wash plugs 2X more with TE buffer
- include BSA (0.1 mg/ml) in enzyme master mix
  - only use molecular grade BSA
  - make aliquots to reduce freeze/thaw cycles
- increase units of enzyme
- use concentrated (40 U/ $\mu$ l vs. 10 U/ $\mu$ l) enzyme to decrease the amount of glycerol in the master mix




➤ **Do not mark during analysis**

# Troubleshooting “ghost” bands

- **Possible solutions, continued**


- confirm water bath temperature is correct
- follow protocol closely – varies for organism and/or enzyme
  - use appropriate buffer
  - use suggested units of enzyme
  - use suggested incubation time
  - use suggested incubation temp
- new vial and/or lot and/or vendor of enzyme and/or buffer
- avoid high-fidelity and fast digest enzymes
- use different spatula for plug slices from each enzyme

H9812 double digest – enzyme not inactivated prior to loading gel



# Troubleshooting “ghost” bands – reagent quality



- Intermittent issues with ghost bands, incomplete digestion
- 10X TBE made in-house
- Change in production of EDTA from supplier
  - same order number but no longer labeled as “*For Molecular Use*”
- Took several months to identify the issue
- Keep track of date, lots, etc... of all reagents



# Effect of Water Quality



- **Difficult to replicate issues due to water quality**
- **Make sure the quality of water used for reagents is the highest possible quality**
  - change filters in water system on regular basis
  - **avoid** storing water in plastic carboys
- **Use **sterile ultrapure** (Type I or Reagent Grade) water for all master mixes**
- ****Non-sterile ultrapure** water can be used for:**
  - buffer used to make gels
  - electrophoresis running buffer
  - de-staining gels
- **Never use tap water**

# Effect of Water Quality



- Sudden, dramatic, but intermittent failures; would seem to resolve and then reappear
- Tested and eliminated all possibilities – organism, enzyme, buffers, etc...
- Water system was cleaned, replaced and gels have looked beautiful since then

# Foreign Particles or Bubbles Embedded in Plugs



# Common Questions

## ❑ **Slanted gels, run length, normalization**

- Power supply – building, dedicated line
- Pump flow rate – check rate annually (at least),  $\sim 70 \approx 1\text{L} / \text{minute}$ , if less then adjust setting, reverse flow to flush out agarose pieces, air bubbles or “kinks” in tubing
- Chiller – temperature accuracy (probe can go bad)
- TBE – variations between vendors and in-house, volume  $\sim 2.2\text{L}$
- Water – use the highest grade possible, not sterile unless indicated
  - Milli-Q is great, DI or “polished” water may be OK if good quality going in
  - quality may change over time – new building, aging filter, storage in plastic carboys may leach chemicals, etc...
- Humidity, air temp, ventilation, neighboring equipment
- Agarose – different lot, different kelp, dirty glassware



# Normalization

✓ **Always check normalization – if bad, do not upload**

❑ **Programming error**

- Accidentally programmed 30kb – 600kb, instead of 700kb, etc...
- *Salmonella* isolates run with *E. coli* conditions and vice versa cannot be analyzed – must re-run with correct conditions

❑ **Wrong agarose**


- “Megabase (Bio-Rad lot# 45100031) does **not** separate large molecular weight bands correctly and **cannot** be used
- **Only** SeaKem Gold (Lonza) can be used for gel running agarose

❑ **Instrument malfunction**

- Electrophoresis box or “brain” may go bad
- Swap out boxes and brains to isolate the problem
- Contact Bio-Rad to have equipment repaired



# Bad Normalization – wrong conditions

- ❑ Gel was run with “*Salmonella*” conditions and analyzed in *Salmonella* database
  - Dark distortion bars indicate poor normalization
  - Re-run with 2.16 s – 63.8 s (30 kb – 700 kb)




# Bad Normalization – wrong conditions

- ❑ Gel was run with “*Salmonella*” conditions and analyzed in *E. coli* database
  - Light distortion bars indicate proper normalization
  - Strongly suggests 2.16 s – 54.17 s (*E. coli* O157) conditions




# Factors Causing Lane Curvature


- ❑ Gel not level when poured
- ❑ Electrophoresis chamber was not level
- ❑ Gel was not flush against black platform
- ❑ Buffer not flowing evenly or not enough buffer in electrophoresis chamber
  - CHEF Mapper requires ~2.2 L in chamber, older models require ~2.0 L
  - clean and check for agarose particles
- ❑ Equipment putting off heat near lines
- ❑ Temperature of buffer fluctuated during run, more than  $14^{\circ}\text{C} \pm 2^{\circ}\text{C}$
- ❑ Broken electrode(s)
- ❑ Mapper “personality”
- ❑ Electrophoresis chamber needs to be serviced









dramatic



subtle



# Unique, sporadic, or chronic issue?



**What to consider** → electrical supply, power interruption, dedicated line, power surge, programming mistake, faulty instrument



**Send Your Troubleshooting  
Questions via e-mail to:**

**pfge@cdc.gov**

**Write “Troubleshooting” in the  
subject line of e-mail message**